Bacillol® 30 Tissues

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Recommendations on spray-wipe disinfection

Inanimate surfaces have been shown to be a major source of contamination and infection (1). To eliminate surfaces that are small and difficult to access as source of contamination and infection, the wipe disinfection is ideally complemented by a spray-wipe disinfection or spray disinfection (2, 3, 4) with an alcohol-based rapid disinfectant. In doing so, it is imperative to apply the products correctly:

- Always prefer a wipe disinfection over the spray or spray-wipe disinfection, as it prevents the formation of aerosols and ensures best possible wetting
- When spraying, wipe afterwards, if possible, to ensure complete wetting (spray-wipe disinfection)
- Limit the simple spray disinfection to areas that cannot be disinfected by using the wipe or spray-wipe procedure
- To largely exclude the risk of inhaling aerosols during the spray disinfection, apply product directly to a dry cloth or spray surfaces from a short distance

- If possible, use alcohol-based products that contain no additives, as these dry without leaving any residue
- Another alternative are alcohol-based disinfection foams that do not form any aerosols during spraying

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HARTMANN-ScandiCare AB